
AERONEBULIZER

family 

MANUALE PER L'UTENTE - USER MANUAL - MANUAL DEL USUARIO - BEDIENUNGSANLEITUNG

COMPONENTI

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Globalcare Medical Technology Co., Ltd.
7th Building, 39 Middle Industrial Main Road, European Industrial Zone, Xiaolan Town, 528415 Zhongshan City, Guangdong Province, REPUBBLICA POPOLARE CINESE
Telefono: +86 760 22589901

Donawa Lifescience
Piazza Albania, 10, 00153 Roma Italia
<https://www.donawa.com/wli/main/contatti.index>

PRODECO PHARMA Srl Unipersonale
Via G.Toniolo, 25/1 - 31033
Castelfranco Veneto (TV) - Italy
www.prodecopharma.com
box@prodecopharma.com

GCE813-IM-MDR-2023-07-10-Rev00

Prodeco Pharma Rev.00

**IMPORTANT:**

Read the instructions carefully before use.

Intended purpose:

The device is an aerosoltherapy system suitable for domestic use. The device is designed for the production of compressed air to operate a small volume nebulizer for the production of medical aerosol for lung respiratory disorders.

Patient population:

The device is intended for use with children from 2 years old, adolescent and adult patients.

Intended users:

The user is an adult person without specific knowledge or professional ability. The patient is the intended operator except in case of children and patients that requires special assistance.

Indications:

Lung acute or chronic diseases of respiratory disorders organs, or inflammation of the upper respiratory system.

Contraindications:

The device is not indicated to use with quick-relief medications during life-threatening asthma attack. No contraindications exist to the administration of aerosols by inhalation. Contraindications related to the medicine used must be checked on the medicine package leaflet. Consult your physician in case of doubts.

**WARNINGS:**

- Use the device only as described in this manual and therefore as an aerosoltherapy system, following the indications of your doctor. Any use different from the intended one is to be considered improper and hence dangerous; the manufacturer cannot be held liable for any damage caused by improper, incorrect and/or unreasonable use, or if the equipment is connected to electrical installations which do not comply with current safety regulations.
- Keep the user manual for future reference.
- Do not operate the unit in presence of any anaesthetic mixture inflammable with oxygen or nitrogen protoxide.
- The correct functioning of the equipment can be affected by electromagnetic interferences which exceed the limits indicated by the European standards in force. In case this device interferes with other electrical devices, move it and plug it to a different power socket.
- In case of failure and/or malfunction, read the "POSSIBLE PROBLEMS AND HOW TO SOLVE THEM" in the user manual. Do not handle or open the compressor housing.
- For repair operations address only to a technical service centre authorized by the manufacturer and require the use of original spare parts. The non-observation of the above mentioned indications can compromise the device safety.
- Comply with the safety regulations concerning the electrical devices and in particular:
 - use only original accessories and components;
 - never submerge the unit in water;
 - never wet the device, it is not protected against water penetration;
 - never touch the unit with wet or moist hands;
 - do not leave the unit exposed to the weather elements;
 - place the unit on a stable and horizontal surface during its operation;
 - the use of this device by children and disabled requires always the close supervision by an adult with full mental faculties;
 - do not pull the power cord or the device itself to unplug it from the power socket;
 - the power plug is the separation element from the grid power; keep the plug accessible when the device is in use.
- Before plugging in the device, make sure that the electrical rating, shown on the rating plate on the bottom of the unit, corresponds to the mains rating.
- In case the power plug provided with the device does not fit your wall socket, address to qualified personnel for the replacement of the plug with a suitable one. In general, the use of adapters, simple or multiple, and/or extension cables is not recommended. If their use is indispensable, it is necessary to use types complying with safety regulations, paying however attention that they do not exceed the maximum power limits, indicated on adapters and extension cables.
- Do not leave the unit plugged in when not in use; unplug the device from the wall socket when it is not operated.

- The installation must be carried out according to the instructions of the manufacturer. An improper installation can cause damage to persons, animals or things, for which the manufacturer cannot be held responsible.
- The power cord of this device cannot be replaced by the user. In case of a power cord damage, address to a technical service center authorized by the manufacturer for its replacement.
- The power supply cord should always be fully unwound in order to prevent dangerous overheating.
- Before performing any maintenance or cleaning operation, turn off the device and disconnect the plug from the main supply.
- Some parts of the unit are so small that they may be swallowed by children; keep the equipment out from children's reach.
- If you decide not to use the device any longer, it is recommended to dispose of it according to the current regulations.
- Make sure to:
 - use this device only with medicines prescribed by your doctor;
 - make the treatment using only the accessory recommended by your doctor depending on the pathology;
 - use the nosepiece accessory only if expressly indicated by your doctor and paying attention to NEVER introduce the bifurcations in the nose, but only bringing them as close as possible.
- Check in the medicine package leaflet for possible contraindications for use with common aerosoltherapy systems.
- To avoid strangulation and entanglement, keep cable and air tubes out of reach of young children.
- Do not position the equipment so that it is difficult to operate the disconnection device.
- Nebulizer and accessories are single patient reuse. Device is multi-patient reuse.
- Nebulizing system is not suitable for use in anaesthetic breathing system or a ventilator breathing system. The patient cannot undergo MRI scan while using this device.
- Device is designed to nebulize solution and suspension liquids.
- Do not use the device, if you spot damage or you notice something unusual.
- This device consists of delicate components and must be handled with care. Observe the storage and operating conditions in the chapter "Technical Data".
- The use of this device is not a replacement for visiting the doctor. In case of no health improvement after treatment consult again your doctor.

HOW TO USE THE DEVICE

The device must be checked before each use in order to identify any malfunctions and/or damage caused during transportation and/or storage. During inhalation, sit upright and relaxed to avoid compressing the airways and impairing treatment effectiveness. The accessories must be used only with a single patient, it is not recommended to use them with several patients.

After unpacking the device, check it for visible damages or defects; pay particular attention to cracks in the plastic housing, which may expose electrical components. Check if the accessories are intact.

Before using the device, proceed with the cleaning operations as described in the chapter «CLEANING AND DISINFECTION».

1. Assemble the nebulizer kit . Ensure that all parts are complete.
2. Fill the nebulizer with the inhalation solution as per your doctor's instructions. Ensure that you do not exceed the maximum level.
3. Connect the nebulizer with the air tube (11) to the compressor and plug the power cord (2)
4. To start the treatment, set ON/OFF switch (1) into the «I» position.
 - The mouthpiece (10) gives you a better drug delivery to the lungs.
 - Choose between adult (9) or child face mask (8) and make sure that it encloses the mouth and nose area completely.
 - Use all accessories including the nose piece (12) as prescribed by your doctor.
5. During inhalation, sit upright and relaxed at a table and not in an armchair, in order to avoid compressing your respiratory airways and impairing the treatment effectiveness. Do not lie down while inhaling. Stop inhalation if you feel unwell.
6. After completing the inhalation period recommended by your doctor, switch the ON/OFF switch (1) to position «O» to turn off the device and unplug it from the socket.
7. Empty the remaining medication from the nebulizer and clean the device as described in the section «CLEANING AND DISINFECTION ».

**CAUTION!**

The device is designed for intermittent use (30 min. on - 30 min. off) Switch off the device after 30 min. use and wait for another 30 min. before you resume treatment.

1. Inhale the aerosol solution using the prescribed accessory.
2. Once you have finished the treatment, switch the appliance off by pressing the I/O switch
3. Clean the nebulizer and its accessories following the steps in the «CLEANING AND MAINTENANCE» section.
4. The device requires no calibration.
5. Tampering with the device is absolutely forbidden.
6. No modification to the device is permitted.

CLEANING AND DISINFECTION

| Device(s): Aerosoltherapy nebulizer accessories | |
|---|--|
| WARNINGS | <p>Before the first use and after each treatment Follow carefully the cleaning and disinfection instructions of the accessories as they are very important for the device performances and the therapy success. Use original accessories only. Do not clean or disinfect the air tube. The accessories can not be cleaned and disinfected by automated method. Do not boil nor autoclave the masks.</p> |
| Limitations on processing | <p>The nebulizer must be replaced after a long period of inactivity, in case it shows deformations or breakings, or when the nebulizer nozzle is obstructed by dry medicine, dust, etc. We recommend replacing the nebulizer after a period between 6 months and 1 year, depending on the usage. The maximum times of cleaning and disinfection of the nebulizer are 360 times. We recommend replacing the nasal wash once a 1 year. The maximum times of cleaning and disinfection of the nasal wash are 360 times. The maximum times of cleaning and disinfection of the mask, the nosepiece and the mouthpiece are 360 times.</p> |
| INSTRUCTIONS | |
| Preparation before cleaning | <p>Detach the air tube from the nebulizer. Detach the mouthpiece or nosepiece or mask from the nebulizer. Make sure that all residual volume is removed from the nebulizer. Disassemble the nebulizer by turning counterclockwise the top and remove the medicine conduction cone. When use the nasal wash: Detach the nasal wash from the air tube. Disassemble the nasal wash. Note: To avoid microbial growth and residual drug drying, clean and disinfect the accessories immediately after each use.</p> |
| Cleaning: | <p>Briefly rinse all parts in running tap water for at least 10 seconds beforehand. Mix little dish detergent and warm tap water i.e. FAIRY brand dish washing liquid for hand washing in a ratio of 2ml: 1L in a clean container. Immerse the components of disassembled nebulizer and disassembled nasal wash (when use the nasal wash), the mask, the mouthpiece and the nosepiece into the mixed water for about 5 min. And then scrub all the surface of all the components with clean and small brush for at least 8 times. Afterwards rinse all parts thoroughly in running tap water at least 30 seconds, to fully remove any possible dish washing liquid residue.</p> |
| Disinfection | <p>After cleaning, disinfect all of the disassembled parts (only parts that have been cleaned can be disinfected effectively).</p> <p>Boil the disassembled nebulizer and disassembled nasal wash (when use the nasal wash), the mouthpiece and the nosepiece for 5 minutes in boiling tap water.</p> <p>Immerse the mask in a 2% sodium hypochlorite (NaOCl) solution for 15 minutes, (i.e. a solution made from the disinfecting agent Amuchina® or a 2% sodium hypochlorite solution mixed by your pharmacist). Afterwards immerse the mask with sterile water for 3 min, and then rinse all the surface of the mask with other clean sterile water twice, to fully remove any possible residue of the disinfectant solution.</p> |

| INSTRUCTIONS | |
|-----------------------|---|
| Drying | <p>Reassemble the nebulizer components and the nasal wash components (when use the nasal wash), and connect the nebulizer and the nasal wash separately to the air-outlet, switch the device on and let it work for 10-15 minutes.</p> <p>Let all parts dry completely before reassembling and reuse to avoid risk of germ grow. Care needs to be taken not to contaminate the parts after they have been cleaned and disinfected. Contamination can be avoided by good hand washing and not touching the inside sections of the device when laying them out to dry or when reassembling.</p> |
| Inspection | Inspect all product components after each cleaning and disinfection. Replace any broken, misshapen or seriously discolored parts. |
| Packaging | Pack dry parts in a clean and sealed container when not in use. Do NOT pack wet or damp parts. |
| Storage | <p>Storage conditions refer to "TECHNICAL SPECIFICATIONS"</p> <p>Note: re-clean and re-disinfect the parts if they are stored more than one day.</p> |
| Transportation | <p>After cleaning and disinfection, always transport the parts with clean and sealed container.</p> <p>Contamination can be avoided by good hand washing and not touching the inside sections of the parts when take out and re-assemble the parts for use.</p> |

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

**FILTER REPLACEMENT**

In normal usage conditions, the air filter must be replaced approximately after 200 working hours or once a year. We recommends to periodically check the air filter(10-12 treatments) and, if the filter shows a grey or brown colour or is wet, replace it. Remove the filter and replace it with a new one.
Do not attempt to clean the filter to reuse it. The air filter shall not be serviced or maintained while in use with a patient.


Only use original filters. Do not use the device without the filter.

Malfunctions and Actions to take:

The device cannot be switched on

- Ensure the power cord (2) is correctly plugged into the socket.
- Ensure the ON/OFF switch (1) is in the position «I».
- Make sure that the device has been operating within operating limits indicated in this manual (30 min On / 30 min Off). The nebulizer functions poorly or not at all
- Ensure the air tube (11) is correctly connected at both ends.
- Ensure the air tube (11) is not squashed, bent, dirty or blocked. If necessary, replace with a new one.
- Ensure the nebulizer (7) is fully assembled and the colored nebulizer nozzle **7-a** is placed correctly and not obstructed.
- Ensure the required medication has been added.

In the event that the equipment does not resume correctly, please contact **CUSTOMER SERVICE**.

**MAINTENANCE AND REPAIRS**

In case of failure, contact **qualified personnel authorised** by Globalcare. **Never open the device**. The device has no user-serviceable parts within and does not need internal maintenance or lubrication.

NASAL WASH INSTRUCTION

GENERAL INFORMATION

Nasal Wash is a medical device suitable for cleaning nasal cavities, it helps to moisturise the nasal mucosa and creates a micronized jet for the treatment of respiratory problems of the upper airways.

We recommend using the nasal wash to relieve the symptoms of colds, dry nose (nasal dryness), dust exposure, pollen allergies, chronic nasal inflammations or as a result of your doctor's advice.



Only use with isotonic saline solution.

The maximum volume of nasal wash solution container is 8 ml. If in doubt, consult your pharmacist to recommend a correct dilution of the isotonic solution for use with nasal wash.

PRODUCT COMPONENTS

- A) Nasal wash cover
- B) Solution container
- C) Waste disposal container
- D) Tubing connector
- E) Aperture for activating aerosol mist

HOW TO USE THE PRODUCT

1. To remove the cover turn 90° (counter clockwise)
2. Fill container with isotonic solution
3. Close the device cover tightly (clockwise)
4. Connect the device to the main unit with the tube
5. Position nose piece on nostril and inhale through nose
6. Turn on main unit then place your finger over the aperture to start treatment

During use you should breathe in and out slowly through your nose and tilt your head slightly to the opposite side of the affected nostril to allow the saline mist to flow deep into your nostril. You can pause or stop the wash by removing your finger from the aperture.

WARNINGS

- Do not use with drugs, balsamic oils and/or other medicinal substances.
- Do not use with ultrasonic nebulizers

| Item | specification | |
|---|---------------|-----|
| Max Pressure (acceptable range) | 200-320 | kPa |
| Max Pressure (average value) | 220 | kPa |
| Operating Pressure (with calibrated nozzle, 0,50) | 80 - 155 | kPa |
| Protection type IP | IP 21 | - |
| Sound noise level (Distance: 1m, Front direction) max level | ≤58 | dBA |

| Item | specification | |
|--|-----------------|--------|
| Sound noise level - average level | 52 | dBA |
| Equipment conditions of use | 30 ON / 30 OFF | Min |
| Dimensions | 160x140x110 | mm |
| Weight | 1,3±0,2kg | kg |
| Service life | 1000 | ore |
| Power supply | 230V CA 50Hz | -- |
| Current ABSORPTION | ≤ 1000 | mA |
| Neb Rate (average value) | 0,35 | ml/min |
| Neb Rate (acceptable range) | 0,23~0,60 | ml/min |
| Nebulizer kit | MedelJet Clenny | |
| AEROSOL PERFORMANCES ACCORDING TO EN ISO 27427:2019 based on adult ventilatory pattern with Sodium Fluoride (NaF): | | |
| Aerosol output | 0,11 | ml |
| Aerosol output rate | 0,07 | ml/min |
| Percentage of fill volume emitted per min | 3,5 | % |
| Residual Volume | 1,68 | ml |
| Particle size (MMAD) | 3,07 | µm |
| GSD (geometric standard deviation) | 2,03 | -- |
| RF (respirable fraction <5µm) | 75 | % |
| Large particle range (>5µm) | 25 | % |
| Middle particle range (2 to 5µm) | 47,8 | % |
| Small particle range (<3µm) | 27,2 | % |
| Min & Max filling volume: | 2 - 8 | ml |
| Operating Temperature | 10 ~ 40 | °C |
| Operating Humidity | 10 ~ 95 | % RH |
| Operating Pressure | 700 ~ 1060 | hPa |
| Storage Temperature | -25 ~ 70 | °C |
| Storage Humidity | 10 ~ 95 | % RH |
| Storage Pressure | 700 ~ 1060 | hPa |
| Protection against electrical shock | 2 | / |

- Class II device referred to protection against electric shock.
- Nebulizer, mouthpiece, masks, and nosepiece are type BF applied parts

Note: Technical specifications subject to change without notice.

AEROSOL CHARACTERISTICS IN ACCORDANCE WITH REGULATION EN ISO27427

Standards applied:

Electric safety standards CEI EN 60601-1 Electromagnetic compatibility according to CEI EN 60601-1-2 The device is a Class IIa medical device. **Device Complies with European Regulation on Medical Devices EU MDR 2017/745**

SAFE DISPOSAL OF THE DEVICE

| Item | Component | Disposal Description | Photo |
|------|---|--|--|
| 1 | Device | The component mainly includes plastic and electronic components. All comply with RoHS and REACH, and all could be safely disposed. This product is subject to European Directive 2012/19/EU on waste electrical and electronic equipment and is marked accordingly. Never dispose of electronic devices with household waste. Please seek out information about the local regulations with regard to the correct disposal of electrical and electronic products. |  |
| 2 | Accessories (air tube, nebulizer, nosepiece, mouthpiece, masks, nasal wash) | The components are plastic. All comply with ROHS and REACH. Before disposal, all accessories should be cleaned according to manual instructions in the chapter «CLEANING AND DISINFECTION». |  |

DEVICE SYMBOLS

| | | | |
|--|--|---|--|
|  | Catalogue Number |  | Distributor |
|  | ON/OFF |  | Authorized representative |
|  | Alternating current |  | Manufacturer |
|  | Class II |  | Serial Number (LOT/serial number) |
|  | Type BF |  | Date of manufacturing (year, month and day) |
|  | Device Complies with European Regulation on Medical Devices EU MDR 2017/745 |  | Caution |
|  | Medical device |  | Humidity limitation |
|  | Protection against solid foreign objects and against harmful effects due to the ingress of water |  | Temperature limitation |
|  | Follow instruction for use |  | Ambient pressure limitation |
|  | Importer |  | Single patient multiple use (for accessories only) |
|  | Unique device Identifier | | |

WEEE- Waste Electrical and Electronic Equipment

This product is subject to European Directive 2012/19/EU on waste electrical and electronic equipment and is marked accordingly. Never dispose of electronic devices with household waste. Please seek out information about the local regulations with regard to the correct disposal of electrical and electronic products. Correct disposal helps to protect the environment and human health.

POSSIBLE PROBLEMS AND HOW TO SOLVE THEM

Contact an authorized Customer Service Centre in case you need help about the use and the maintenance of the equipment.

PROBLEMS:

- The device does not switch on
- The device does not nebulize or nebulizes weakly

SOLUTIONS:

- Make sure the power plug is firmly fitted to the wall socket.
 - Make sure that the device has been operating within operating limits indicated in this manual (30 min ON / 30 min OFF).
 - Make sure that the ends of the air tube are tightly to the main unit and nebulizer.
 - Verify that the nebulizer is not empty or has been filled with the proper quantity of medicine (max 8 ml). Verify that the nebulizer nozzle is not obstructed.
- Should the device not start working properly again, address to customer service of authorized retailer.

MAINTENANCE AND REPAIRS

In case of failure, address to qualified personnel authorized by manufacturer. Do not open the device in any case. The unit has no user-serviceable parts within and does not need internal maintenance or lubrication.
Please report any serious incident that has occurred in relation to the device injury or adverse event to the local competent authority and to the Manufacturer or to the European Authorised Representative (EC REP)
Vigilance contact points: <https://ec.europa.eu/growth/sectors/medical-devices/contacts/>

Electromagnetic compatibility

Electromagnetic emissions

Guidance and manufacturer's declaration - electromagnetic emissions for all EQUIPMENT and SYSTEMS. The Nebulizer is intended for use in the electromagnetic environment specified below. The customer of the user of the Nebulizer should assure that it is used in such an environment.

| Emission test | Compliance | Electromagnetic environment- guidance |
|------------------------------|------------|--|
| RF emissions CISPR 11 | Group 1 | The Nebulizer uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emission CISPR 11 | Class B | The Nebulizer is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |

Electromagnetic immunity

Guidance and manufacturer's declaration - electromagnetic immunity - for all EQUIPMENT and SYSTEMS. The Nebulizer is intended for use in the electromagnetic environment specified below. The customer of the user of the Nebulizer should assure that it is used in such an environment.

| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidance |
|--|--|-----------------------------|--|
| Electrostatic discharge (ESD) IEC 61000-4-2 | ±8 kV contact ±15 kV air | ±8 kV contact ±15 kV air | Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%. |
| Power frequency (50Hz) magnetic Field IEC 61000-4-8 | 30 A/m 50Hz/60Hz | 30 A/m 50Hz/60Hz | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |
| Proximity magnetic fields (IEC 61000-4-39) | 30kHz, CW, 8A/m 134,2kHz, PM 2,1kHz, 65A/m 13,56MHz, PM 50kHz, 7,5A/m | N/D | The EUT does not contain magnetically sensitive components or circuitry, so this test does not need to be evaluated. |

Guidance and manufacturer's declaration - electromagnetic immunity for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

The Nebulizer is intended for use in the electromagnetic environment specified below. The customer or the user of Nebulizer should assure that it is used in such an environment.

| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidance |
|--|--------------------------|------------------|---|
| Conducted RF IEC 61000-4-6 | 150KHz to 80MHz 3Vrms | N/D | Portable and mobile RF communications equipment should be used no closer to any part of the Nebulizer, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d=1.2\sqrt{P}$ $d=1.2\sqrt{P}$ |
| Radiated RF IEC 61000-4-3 | 10 V/m | 10 V/m | $d=1.2\sqrt{P}$ 80 MHz to 800 GHz $d=2.3\sqrt{P}$ 800 MHz to 2.7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey a should be less than the compliance level in each frequency range b. Interference may occur in the vicinity of equipment marked with the following symbol:  |
| NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. | | | |
| (a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Nebulizer is used exceeds the applicable RF compliance level above, the Nebulizer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Nebulizer. (b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m. | | | |

Recommended separation distances

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM - for EQUIPMENT or SYSTEM that are not LIFESUPPORTING.

Recommended separation distances between portable and mobile RF communications equipment and the Nebulizer. The Nebulizer is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Nebulizer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Nebulizer as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output Power of transmitter (W) | Separation distance according to frequency of transmitter (m) | | | |
|---|---|--|--------------------------------------|---------------------------------------|
| | 150 kHz to 80 MHz(out ISM and amateur radio bands) $d=1.2\sqrt{P}$ | 150 kHz to 80 MHz(in ISM and amateur radio bands) $d=2\sqrt{P}$ | 80 MHz to 800 MHz $d=1.2\sqrt{P}$ | 800 MHz to 2.7 GHz $d=2.3\sqrt{P}$ |
| 0,01 | 0,12 | 0,2 | 0,12 | 0,23 |
| 0,1 | 0,38 | 0,632 | 0,38 | 0,73 |
| 1 | 1,2 | 2 | 1,2 | 2,3 |
| 10 | 3,8 | 6,32 | 3,8 | 7,3 |
| 100 | 12 | 20 | 12 | 23 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Globalcare Medical Technology Co., Ltd.
7th Building, 39 Middle Industrial Main Road, European Industrial Zone, Xiaolan Town, 528415 Zhongshan City, Guangdong Province, PEOPLE'S REPUBLIC OF CHINA.
Phone: +86 760 22589901

Donawa Lifescience
Piazza Albania, 10, 00153 Rome Italy
<https://www.donawa.com/wli/main/contatti.index>

PRODECO PHARMA Srl Unipersonale
Via G.Toniolo, 25/1 - 31033
Castelfranco Veneto (TV) - Italy
www.prodecopharma.com
box@prodecopharma.com

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SÍMBOLOS DEL DISPOSITIVO

| | | | |
|--|--|---|--|
|  | Catalogue Number |  | Distributor |
|  | ON/OFF |  | Authorized representative |
|  | Alternating current |  | Manufacturer |
|  | Class II |  | Serial Number (LOT/serial number) |
|  | Type BF |  | Date of manufacturing (year, month and day) |
|  | Device Complys with European Regulation on Medical Devices EU MDR 2017/745 |  | Caution |
|  | Medical device |  | Humidity limitation |
|  | Protection against solid foreign objects and against harmful effects due to the ingress of water |  | Temperature limitation |
|  | Follow instruction for use |  | Ambient pressure limitation |
|  | Importer |  | Single patient multiple use (for accessotied only) |
|  | Unique device Identifier | | |

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GERÄTESYMBOLS

- Catalogue Number
- ON/OFF
- Alternating current
- Class II
- Type BF
- Device Complys with European Regulation on Medical Devices EU MDR 2017/745
- Medical device
- Protection against solid foreign objects and against harmful effects due to the ingress of water
- Follow instruction for use
- Importer
- Unique device Identifier
- Distributor
- Authorized representative
- Manufacturer
- Serial Number (LOT/serial number)
- Date of manufacturing (year, month and day)
- Caution
- Humidity limitation
- Temperature limitation
- Ambient pressure limitation
- Single patient multiple use (for accessotied only)

WEEE (Waste Electrical and Electronic Equipment; Elektro- und Elektronikgeräte-Abfall)

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Globalcare Medical Technology Co. Ltd.